



VideoTest-Album

Image Database

Windows 9x/NT compatible software

VideoTest- Album Image Database is designed to store images and data (histograms, tables, and text) and perform their quick search. It also provides an ability to generate reports and print them.

VideoTest- Album main functions:

- Import the images from Videotest software, open them from files and insert via clipboard
- Adjustable image compression
- Enter the necessary image information
- User-defined view of the record
- Ability to manage the unlimited number of databases
- Quick search for data of interest
- Ability to create queries
- Ability to generate reports and print them
- Image export to VideoTest software
- Export the images and data information to other Windows applications.

VideoTest- Album has three modes:

Card mode is used to fill new information in. The Card contains an image (images) and data (images, tables, histograms, graphs, text, etc). The user can easily configure the Card view.

Gallery mode displays all data base images. The image size and comment can be adjusted by the user.


Report generator mode is used to create standard reports. The report contains the user-selected fields of the record.



Card



Gallery



Report



P.O. Box 356, St.-Petersburg, 190000 **RUSSIA**
Phone: +7(812) 314-8100, 314-8445
Fax +7(812) 325-6494
e-mail: info@videotest.ru
<http://www.videotest.ru>

Distributed by: