



# VideoTest-Colony Count



## Counting bacteria colonies

**VideoTest-Colony Count** is designed to acquire, separate, count and measure one, two, three or four bacteria colonies. The colony images are captured with color or b/w CCD camera of high resolution. The measure data are processed statistically. The images and data can be saved in VideoTest-Album image database.



### VideoTest-Colony Count main functions:

- Acquisition of the colonia images with CCD camera, scanner, open from disk or load from VideoTest-Album Image Database,
- Image enhancement and annotation,
- Ability to use four pre-installed methodics designed to count and size the bacteria colonies of one, two, three or four types,
- Automatic separation (thresholding) of the bacteria colonies,
- Easy editing the separated colonies,
- Count and measure the bacteria colonies with data display in the spreadsheet,
- Statistical treatment of the measure results,
- Ability to add new methodics of analysis,
- Ability to copy and modify the pre-installed methodics,
- Image and data export to VideoTest-Album 3.0 image database (not included).



P.O. Box 356, St.-Petersburg, 190000 **RUSSIA**  
Phone: +7(812) 314-8100, 314-8445  
Fax +7(812) 325-6494  
e-mail: [info@videotest.ru](mailto:info@videotest.ru)  
<http://www.videotest.ru>

Distributed by: