



Forensic Pro

Finger Print Analysis Image Comparison

Forensic Pro software is designed for capturing and comparing the images (finger prints, tire tracks, shoe prints, bite marks, etc.). Forensic Pro provides the abilities of image measurements, annotations, enhancement, has user-friendly flexible interface and easy-to-learn.

Main functions of **Forensic Pro**:

- B/w or color Image acquisition with CCD / digital cameras (including 10, 12, 16 bit), scanner, open from files, paste from the clipboard,
- Ability to work with series of images and analysis results within one document. Convenient switching between images and data. Save series of images and data to the same document,
- Image enhancement with filters, FFT transformations,
- Image annotation (text, graphics), move, copy, rotate the outlined image area, zoom in and out, graphics exchange between images,
- Overlay images with traces,
- Manual measurements (linear, angle, radiuses, distance between two parallel lines, etc.), count objects. All measure data are displayed in the spreadsheets,
- Automatic measurement of the areas (objects) manually outlined,
- Statistical treatment of the measurement results, data export to MS Excel,
- Numerous abilities to compare images, synchronous image scrolling and tracing overlay,
- "Audit trails" with an ability of recording all operations performed under the images,
- Ability to create the multi-page reports,
- Save and print images and data,
- VideoTesT-Album 4.0 image database included.



VideoTest

P.O. Box 356, St.-Petersburg, 190000 **RUSSIA**
Phone: +7(812) 314-8100, 314-8445
Fax +7(812) 325-6494
E-mail: info@videotest.ru
<http://www.videotest.ru>